

LABMATE PRO VARIABLE VOLUME PIPETTE SERIES

SINGLE CHANNEL PIPETTES							MULTICHANNEL PIPETTES						
MODEL	CAT. NO.	COLOUR	VOLUME (µL)	A (%)	P (%)	FIT TO TIPS	MODEL	CAT. NO.	COLOUR	VOLUME (µL)	A (%)	P (%)	FIT TO TIPS
LMP2*	5661	RED	0.2	± 12.0	± 6.0	10µL	LMP8-10	6281	RED	MIN 0.5	± 10.0	± 8.0	10µL
			1.0	± 2.7	± 1.3					5	± 4.0	± 2.0	
			MAX 2.0	± 1.5	± 0.7		MAX 10	± 2.0	± 1.2				
LMP10	5662	RED	MIN 0.5	± 4.0	± 2.8		LMP8-50	6282	YELLOW	MIN 5	± 4.0	± 2.5	200µL
			5	± 1.0	± 0.6					25	± 3.0	± 1.2	
MAX 10	± 0.5	± 0.4	MAX 50	± 1.6	± 0.6								
LMP20	5663	YELLOW	MIN 2	± 3.0	± 1.5		LMP8-200	6283	YELLOW	MIN 20	± 3.0	± 1.5	
			10	± 1.0	± 0.5					100	± 1.5	± 0.8	
			MAX 20	± 0.8	± 0.3		MAX 200	± 1.0	± 0.6				
LMP50	5667	YELLOW	MIN 5	± 2.5	± 2.0		LMP8-300	6284	GREEN	MIN 50	± 1.6	± 1.5	300µL
			25	± 1.0	± 0.6	150				± 1.2	± 1.0		
MAX 50	± 0.8	± 0.4	MAX 500	± 1.0	± 0.6								
LMP100	5664	YELLOW	MIN 10	± 1.6	± 0.8								
			50	± 0.8	± 0.24								
MAX 100	± 0.8	± 0.2											
LMP200	5665	YELLOW	MIN 20	± 1.2	± 0.6								
			100	± 0.8	± 0.25								
MAX 200	± 0.6	± 0.2											
LMP250	5660	GREEN	MIN 50	± 1.0	± 0.4								
			125	± 0.8	± 0.3								
MAX 250	± 0.6	± 0.3											
LMP1000	5666	BLUE	MIN 100	± 1.6	± 0.4								
			500	± 0.7	± 0.2								
MAX 1000	± 0.6	± 0.15											
LMP5000	5668	WHITE	MIN 500	± 1.2	± 0.5								
			2500	± 0.6	± 0.2								
MAX 5000	± 0.5	± 0.15											
LMP10000	5669	WHITE	MIN 1000	± 2.5	± 0.6								
			5000	± 0.8	± 0.3								
MAX 10000	± 0.5	± 0.2											

* LMP2 – minimum volume 0.1 µl

© 2015 PZ HTL S.A. All rights reserved.