

MPW4+ Microplate Washer



Features

- 1- Fully automatic
- 2- User friendly operation system with large LCD display
- 3- Washes flat, U and V-bottom strips and plate
- 4- 12-way or 8-way manifold, automatic and manual positioning
- 5- Large memory, up to 100 wash protocols can be programmed
- 6- Multi-channel, 1 or 3 (optional) for wash, 1 for rinse and 1 for waste
- 7- Plate shaking function, time and speed are adjustable
- 8- Two 96-well plate incubation positions (optional)
- 9- Liquid level detection and alert function
- 10- Easy maintenance

System Specifications

Plate Mode	96-/48-well plate or strip
Residual Volume	31 μ l
Washing Volume	10~3000 μ l/well, 1 μ l by step
Washing Cycle	Maximum 99 cycles
Strip Washing	1 to 12 strips programmable
Soaking Time	0~24hours
Shaking Time	0~24hours
Display	Large LCD
Wash Channels	5 channel bottles, 1 or 3 (optional) for wash, 1 for water(rinse) and 1 for waste
Bottles Capacity	2 Liters each bottle
Memory	100 washing protocols
Interface	RS-232
Power Supply	AC 220V \pm 10%, 50HZ
Net Weight	13 KG
Dimension LxWxH(mm)	450x390x190

Hiperion



Medizintechnik GmbH & Co.KG
Carl-Zeiss-Str.18 D-6332 Roedermark,Germany

